

For Research Use Only

Cytokeratin 14 Recombinant antibody, PBS Only



Catalog Number: 83379-1-PBS

Basic Information

Catalog Number:

83379-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0188

GenBank Accession Number:

BC002690

GeneID (NCBI):

3861

UNIPROT ID:

P02533

Full Name:

keratin 14

Calculated MW:

472 aa, 52 kDa

Purification Method:

Protein A purification

CloneNo.:

1H22

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat skin tissue were subjected to SDS PAGE followed by western blot with 83379-1-RR (KRT14 antibody) at dilution of 1:1400 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83379-1-PBS in a different storage buffer formulation.