

For Research Use Only

DSC3 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83386-4-PBS



Basic Information

Catalog Number:

83386-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34648

GenBank Accession Number:

NM_001941

GeneID (NCBI):

1825

UNIPROT ID:

Q14574

Full Name:

desmocollin 3

Calculated MW:

100 kDa

Purification Method:

Protein A purification

CloneNo.:

240292H4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00421-1, DSC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83386-2-PBS. Detection antibody: 83386-4-PBS. Standard: Ag34648. Range: 0.312-40 ng/mL



Cytometric bead array standard curve of MP00421-3, DSC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83386-4-PBS. Detection antibody: 83386-3-PBS. Standard: Ag34648. Range: 0.625-20 ng/mL



Biolayer Interferometry (BLI) kinetic assay of 83386-4-PBS against Human DSC3 was performed. The affinity constant is 3.65 nM.