

For Research Use Only

# C9orf25 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83387-1-PBS

## Basic Information

Catalog Number:

83387-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23771

GenBank Accession Number:

BC041009

GeneID (NCBI):

203259

UNIPROT ID:

Q8IW50

Full Name:

chromosome 9 open reading frame 25

Purification Method:

Protein A purification

CloneNo.:

240151E11

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00426-1, C9orf25 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83387-1-PBS. Detection antibody: 83387-3-PBS. Standard: Ag23771. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP00426-2, C9orf25 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83387-2-PBS. Detection antibody: 83387-1-PBS. Standard: Ag23771. Range: 0.625-40 ng/mL.