

For Research Use Only

TMEM189 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83388-1-PBS

Basic Information

Catalog Number:

83388-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34323

GenBank Accession Number:

BC142966

GeneID (NCBI):

387521

UNIPROT ID:

A5PLL7

Full Name:

transmembrane protein 189

Purification Method:

Protein A purification

CloneNo.:

240376B5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00431-2, TMEM189 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83388-4-PBS. Detection antibody: 83388-1-PBS. Standard: Ag34323. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00431-3, TMEM189 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83388-5-PBS. Detection antibody: 83388-1-PBS. Standard: Ag34323. Range: 0.78-100 ng/mL