

For Research Use Only

SIAH1 Recombinant antibody, PBS Only

Catalog Number: 83389-3-PBS



Basic Information

Catalog Number:

83389-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34618

GenBank Accession Number:

BC035562

GeneID (NCBI):

6477

UNIPROT ID:

Q8IUQ4

Full Name:

seven in absentia homolog 1
(Drosophila)

Calculated MW:

282 aa, 31 kDa

Purification Method:

Protein A purification

CloneNo.:

240377G5

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assay of 83389-3-PBS against Human SIAH1 was performed. The affinity constant is 110 nM.