

For Research Use Only

NCOA4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83394-2-PBS



Basic Information

Catalog Number:

83394-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15277

GenBank Accession Number:

BC012736

GeneID (NCBI):

8031

UNIPROT ID:

Q13772

Full Name:

nuclear receptor coactivator 4

Calculated MW:

70 kDa

Purification Method:

Protein A purification

CloneNo.:

240169F9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00435-1, NCOA4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83394-4-PBS. Detection antibody: 83394-2-PBS. Standard: Ag15277. Range: 3.13-100 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83394-2-PBS against Human NCOA4 was performed. The affinity constant is 0.177 nM.