

For Research Use Only

# FAM105B Recombinant antibody, PBS Only (Detector)

Catalog Number: 83395-1-PBS



## Basic Information

Catalog Number:

83395-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22302

GenBank Accession Number:

BC007706

GeneID (NCBI):

90268

UNIPROT ID:

Q96BN8

Full Name:

family with sequence similarity 105, member B

Calculated MW:

352 aa, 40 kDa

Purification Method:

Protein A purification

CloneNo.:

240267D5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00438-3, FAM105B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83395-2-PBS. Detection antibody: 83395-1-PBS. Standard: Ag22302. Range: 3.13-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83395-1-PBS against Human FAM105B was performed. The affinity constant is 0.39 nM.