

For Research Use Only

FAM105B Recombinant antibody, PBS Only (Detector)

Catalog Number: 83395-3-PBS



Basic Information

Catalog Number:

83395-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22302

GenBank Accession Number:

BC007706

GeneID (NCBI):

90268

UNIPROT ID:

Q96BN8

Full Name:

family with sequence similarity 105, member B

Calculated MW:

352 aa, 40 kDa

Purification Method:

Protein A purification

CloneNo.:

240267E2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00438-1, FAM105B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83395-2-PBS. Detection antibody: 83395-3-PBS. Standard: Ag22302. Range: 3.13-100 ng/mL.