

For Research Use Only

# MPHOSPH9 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83397-2-PBS

## Basic Information

Catalog Number:

83397-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22228

GenBank Accession Number:

BC130440

GeneID (NCBI):

10198

UNIPROT ID:

Q99550

Full Name:

M-phase phosphoprotein 9

Calculated MW:

1031 aa, 116 kDa

Purification Method:

Protein A purification

CloneNo.:

240182F6

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00416-3, MPHOSPH9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83397-2-PBS. Detection antibody: 83397-3-PBS. Standard: Ag22228. Range: 1.25-40 ng/mL