

For Research Use Only

# Fas/CD95 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83400-3-PBS



## Basic Information

Catalog Number:

83400-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_000043.6

GeneID (NCBI):

355

UNIPROT ID:

P25445-1

Full Name:

Fas (TNF receptor superfamily, member 6)

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

240340E8

## Applications

Tested Applications:

Cytometric bead array, Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


# Selected Validation Data



Sandwich ELISA standard curve of MP00430-2, human Fas/CD95 Recombinant Matched Antibody Pair - PBS only. 83400-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0877. 83400-3-PBS was HRP conjugated as the detection antibody. Range: 0.39-25 ng/mL.



Serum of fourteen individual healthy human donors was measured. The Fas/TNFRSF6/CD95 concentration of detected samples was determined to be 2.33 ng/mL with a range of 0.01 - 6.29 ng/mL.



Cytometric bead array standard curve of MP00430-2, Fas/CD95 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83400-2-PBS. Detection antibody: 83400-3-PBS. Standard: Eg0877. Range: 1.25-40 ng/mL.