

For Research Use Only

# FGFR3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83402-1-PBS



## Basic Information

Catalog Number:

83402-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

2261

UNIPROT ID:

P22607-1

Full Name:

fibroblast growth factor receptor 3

Calculated MW:

88kd

Purification Method:

Protein A purification

CloneNo.:

240378A1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00407-1, FGFR3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83402-1-PBS. Detection antibody: 83402-2-PBS. Standard: Eg1040. Range: 3.125-50 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83402-1-PBS against Human FGFR3 was performed. The affinity constant is 0.315 nM.