

For Research Use Only

SREBF2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83403-4-PBS

Basic Information

| | | | | | |
|---------------------------|-------------|---------------------------|--|----------------------|------------------------|
| Catalog Number: | 83403-4-PBS | GenBank Accession Number: | BC056158 | Purification Method: | Protein A purification |
| Size: | 1 mg/ml | GeneID (NCBI): | 6721 | CloneNo.: | 240451C8 |
| Source: | Rabbit | UNIPROT ID: | Q12772 | | |
| Isotype: | IgG | Full Name: | sterol regulatory element binding transcription factor 2 | | |
| Immunogen Catalog Number: | AG28205 | Calculated MW: | 124 kDa | | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00408-3, SREBF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83403-4-PBS. Detection antibody: 83403-1-PBS. Standard: Ag28205. Range: 0.313-40 ng/mL.