

For Research Use Only

SLC25A17 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83404-3-PBS

Basic Information

Catalog Number:

83404-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34310

GenBank Accession Number:

BC012998

GeneID (NCBI):

10478

UNIPROT ID:

O43808

Full Name:

solute carrier family 25
(mitochondrial carrier; peroxisomal
membrane protein, 34kDa), member
17

Calculated MW:

307 aa, 35 kDa

Purification Method:

Protein A purification

CloneNo.:

240382F5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00409-2, SLC25A17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83404-3-PBS. Detection antibody: 83404-2-PBS. Standard: Ag34310. Range: 2.5-80 ng/mL.