

For Research Use Only

# Mouse IL-22 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83405-1-PBS



## Basic Information

Catalog Number:

83405-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0405

GenBank Accession Number:

NM\_016971

GeneID (NCBI):

50929

UNIPROT ID:

Q9JJY9

Full Name:

interleukin 22

Calculated MW:

20kd

Purification Method:

Protein A purification

CloneNo.:

240251A12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00410-2, MOUSE IL-22 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83405-1-PBS. Detection antibody: 83405-2-PBS. Standard: Eg0405. Range: 0.625-40 ng/mL