

For Research Use Only

# SLC2A6 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83406-2-PBS

## Basic Information

Catalog Number:	83406-2-PBS	GenBank Accession Number:	BC013740	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	11182	CloneNo.:	240388F10
Source:	Rabbit	UNIPROT ID:	Q9UGQ3		
Isotype:	IgG	Full Name:	solute carrier family 2 (facilitated glucose transporter), member 6		
Immunogen Catalog Number:	AG18481	Calculated MW:	55 kDa		

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00414-2, SLC2A6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83406-2-PBS. Detection antibody: 83406-3-PBS. Standard: Ag18481. Range: 3.13-100 ng/mL