

For Research Use Only

SLC2A6 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83406-3-PBS

Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
83406-3-PBS	BC013740	Protein A purification
Size:	GeneID (NCBI):	CloneNo.:
1 mg/ml	11182	240388F4
Source:	UNIPROT ID:	
Rabbit	Q9UGQ3	
Isotype:	Full Name:	
IgG	solute carrier family 2 (facilitated glucose transporter), member 6	
Immunogen Catalog Number:	Calculated MW:	
AG18481	55 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00414-1, SLC2A6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83406-1-PBS. Detection antibody: 83406-3-PBS. Standard: Ag18481. Range: 1.56-100 ng/mL



Cytometric bead array standard curve of MP00414-2, SLC2A6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83406-2-PBS. Detection antibody: 83406-3-PBS. Standard: Ag18481. Range: 3.13-100 ng/mL