

For Research Use Only

# Collagen Type XXVIII Recombinant antibody, PBS Only (Capture)

Catalog Number: 83408-1-PBS



## Basic Information

Catalog Number:

83408-1-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21459

GenBank Accession Number:

BC136892

GeneID (NCBI):

340267

UNIPROT ID:

Q2UY09

Full Name:

collagen, type XXVIII, alpha 1

Calculated MW:

1125 aa, 117 kDa

Purification Method:

Protein A purification

CloneNo.:

240387B1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00423-1, COL28A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83408-1-PBS. Detection antibody: 83408-2-PBS. Standard: Ag21459. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP00423-2, COL28A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83408-1-PBS. Detection antibody: 83408-3-PBS. Standard: Ag21459. Range: 0.625-80 ng/mL.