

For Research Use Only

TRNP1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83410-1-PBS

Basic Information

Catalog Number:	83410-1-PBS	GenBank Accession Number:	BC069216	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	388610	CloneNo.:	240331B12
Source:	Rabbit	UNIPROT ID:	Q6NT89		
Isotype:	IgG	Full Name:	TMF1-regulated nuclear protein 1		
Immunogen Catalog Number:	AG22504	Calculated MW:	227 aa, 23 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00432-2, TRNP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83410-2-PBS. Detection antibody: 83410-1-PBS. Standard: Ag22504. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00432-3, TRNP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83410-3-PBS. Detection antibody: 83410-1-PBS. Standard: Ag22504. Range: 1.25-80 ng/mL