

For Research Use Only

C13orf18 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83411-1-PBS



Basic Information

Catalog Number:

83411-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15474

GenBank Accession Number:

BC032311

GeneID (NCBI):

80183

UNIPROT ID:

Q9H714

Full Name:

chromosome 13 open reading frame

18

Calculated MW:

662 aa, 73 kDa

Purification Method:

Protein A purification

CloneNo.:

240259C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00398-3, C13orf18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83411-5-PBS. Detection antibody: 83411-1-PBS. Standard: Ag15474. Range: 0.625-40 ng/mL.