

For Research Use Only

# C13orf18 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83411-5-PBS



## Basic Information

Catalog Number:

83411-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15474

GenBank Accession Number:

BC032311

GeneID (NCBI):

80183

UNIPROT ID:

Q9H714

Full Name:

chromosome 13 open reading frame  
18

Calculated MW:

662 aa, 73 kDa

Purification Method:

Protein A purification

CloneNo.:

240259F12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00398-1, C13orf18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83411-5-PBS. Detection antibody: 83411-3-PBS. Standard: Ag15474. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00398-3, C13orf18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83411-5-PBS. Detection antibody: 83411-1-PBS. Standard: Ag15474. Range: 0.625-40 ng/mL