

For Research Use Only

# SLC6A7 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83412-1-PBS



## Basic Information

Catalog Number:

83412-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34413

GenBank Accession Number:

BC093785

GeneID (NCBI):

6534

UNIPROT ID:

Q99884

Full Name:

solute carrier family 6  
(neurotransmitter transporter, L-  
proline), member 7

Calculated MW:

636 aa, 71 kDa

Purification Method:

Protein A purification

CloneNo.:

240395D4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00399-2, SLC6A7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83412-1-PBS. Detection antibody: 83412-2-PBS. Standard: Ag34413. Range: 1.25-80 ng/mL.