

For Research Use Only

KIAA1715 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83418-1-PBS



Basic Information

Catalog Number:

83418-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34669

GenBank Accession Number:

BC105132

GeneID (NCBI):

80856

UNIPROT ID:

Q9C0E8

Full Name:

KIAA1715

Purification Method:

Protein A purification

CloneNo.:

240211A10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00417-1, KIAA1715 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83418-1-PBS. Detection antibody: 83418-3-PBS. Standard: Ag34669. Range: 1.25-40 ng/mL



Cytometric bead array standard curve of MP00417-2, KIAA1715 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83418-2-PBS. Detection antibody: 83418-1-PBS. Standard: Ag34669. Range: 2.5-80 ng/mL