

For Research Use Only

Mouse CD117/c-Kit Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83419-1-PBS



Basic Information

Catalog Number:

83419-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0883

GenBank Accession Number:

NM_001122733.1

GeneID (NCBI):

16590

UNIPROT ID:

P05532-1

Full Name:

kit oncogene

Calculated MW:

109 kDa

Purification Method:

Protein A purification

CloneNo.:

240246D11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00418-1, MOUSE CD117/c-kit Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83419-1-PBS. Detection antibody: 83419-2-PBS. Standard: Eg0883. Range: 0.781-25 ng/mL



Cytometric bead array standard curve of MP00418-2, MOUSE CD117/c-kit Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83419-3-PBS. Detection antibody: 83419-1-PBS. Standard: Eg0883. Range: 6.25-100 ng/mL