

For Research Use Only

# TMEM126A Recombinant antibody, PBS Only (Detector)



Catalog Number: 83440-1-PBS

## Basic Information

Catalog Number: 83440-1-PBS	GenBank Accession Number: BC007875	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 84233	CloneNo.: 240425A7
Source: Rabbit	UNIPROT ID: Q9H061	
Isotype: IgG	Full Name: transmembrane protein 126A	
Immunogen Catalog Number: AG14705	Calculated MW: 195 aa, 22 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00449-1, TMEM126A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83440-2-PBS. Detection antibody: 83440-1-PBS. Standard: Ag14705. Range: 6.25-100 ng/mL.