

For Research Use Only

SLC2A10 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83441-2-PBS

Basic Information

Catalog Number: 83441-2-PBS	GenBank Accession Number: BC101657	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 81031	CloneNo.: 240335H8
Source: Rabbit	UNIPROT ID: O95528	
Isotype: IgG	Full Name: solute carrier family 2 (facilitated glucose transporter), member 10	
Immunogen Catalog Number: AG14210	Calculated MW: 541 aa, 57 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00450-1, SLC2A10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83441-2-PBS. Detection antibody: 83441-1-PBS. Standard: Ag14210. Range: 0.78-100 ng/mL.