

For Research Use Only

TMEM119 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83448-1-PBS



Basic Information

Catalog Number:

83448-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_181724.2

GeneID (NCBI):

338773

UNIPROT ID:

Q4V9L6

Full Name:

transmembrane protein 119

Calculated MW:

29 kDa

Purification Method:

Protein A purification

CloneNo.:

240433B9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00444-2, TMEM119 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83448-1-PBS. Detection antibody: 83448-2-PBS. Standard: Eg0796. Range: 0.625-40 ng/mL



Cytometric bead array standard curve of MP00444-3, TMEM119 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83448-1-PBS. Detection antibody: 83448-4-PBS. Standard: Eg0796. Range: 1.25-40 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83448-1-PBS against Human TMEM119 was performed. The affinity constant is 9.53 nM.