

For Research Use Only

# PLXDC2 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83449-2-PBS



## Basic Information

Catalog Number:

83449-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2940

GenBank Accession Number:

BC012885

GeneID (NCBI):

84898

UNIPROT ID:

Q6UX71

Full Name:

plexin domain containing 2

Calculated MW:

480 aa, 54 kDa

Purification Method:

Protein A purification

CloneNo.:

240174F1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00445-1, PLXDC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83449-3-PBS. Detection antibody: 83449-2-PBS. Standard: Ag2940. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00445-3, PLXDC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83449-2-PBS. Detection antibody: 83449-1-PBS. Standard: Ag2940. Range: 0.313-40 ng/mL