

For Research Use Only

# ABCG8 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83450-2-PBS

## Basic Information

Catalog Number:	83450-2-PBS	GenBank Accession Number:	BC113657	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	64241	CloneNo.:	240293E4
Source:	Rabbit	UNIPROT ID:	Q9H221		
Isotype:	IgG	Full Name:	ATP-binding cassette, sub-family G (WHITE), member 8		
Immunogen Catalog Number:	AG18526	Calculated MW:	673 aa, 76 kDa		

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00446-2, ABCG8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83450-2-PBS. Detection antibody: 83450-4-PBS. Standard: Ag18526. Range: 1.25-80 ng/mL.