

For Research Use Only

GNB3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83452-3-PBS

Basic Information

Catalog Number:	83452-3-PBS	GenBank Accession Number:	BC000115	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	2784	CloneNo.:	240281C5
Source:	Rabbit	UNIPROT ID:	P16520		
Isotype:	IgG	Full Name:	guanine nucleotide binding protein (G protein), beta polypeptide 3		
Immunogen Catalog Number:	AG0127	Calculated MW:	37 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00443-1, GNB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83452-3-PBS. Detection antibody: 83452-2-PBS. Standard: Ag0127. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00443-2, GNB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83452-3-PBS. Detection antibody: 83452-1-PBS. Standard: Ag0127. Range: 0.78-100 ng/mL