

For Research Use Only

HGF Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83453-2-PBS

Basic Information

Catalog Number:

83453-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25478

GenBank Accession Number:

BC130286

GeneID (NCBI):

3082

UNIPROT ID:

P14210

Full Name:

hepatocyte growth factor
(hepapoietin A; scatter factor)

Purification Method:

Protein A purification

CloneNo.:

240383E11

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00454-1, HGF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83453-2-PBS. Detection antibody: 83453-1-PBS. Standard: Ag25478. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00454-2, HGF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83453-3-PBS. Detection antibody: 83453-2-PBS. Standard: Ag25478. Range: 0.313-40 ng/mL