

For Research Use Only

# IFITM2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83455-1-PBS

## Basic Information

Catalog Number:	83455-1-PBS	GenBank Accession Number:	BC009696	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	10581	CloneNo.:	240318B7
Source:	Rabbit	UNIPROT ID:	Q01629		
Isotype:	IgG	Full Name:	interferon induced transmembrane protein 2 (1-8D)		
Immunogen Catalog Number:	AG3451	Calculated MW:	132 aa, 15 kDa		

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00456-3, IFITM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83455-4-PBS. Detection antibody: 83455-1-PBS. Standard: Ag3451. Range: 0.78-100 ng/mL