

For Research Use Only

PMS1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83456-8-PBS



Basic Information

Catalog Number:

83456-8-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28132

GenBank Accession Number:

BC096331

GeneID (NCBI):

5378

UNIPROT ID:

P54277

Full Name:

PMS1 postmeiotic segregation
increased 1 (S. cerevisiae)

Purification Method:

Protein A purification

CloneNo.:

240408C5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00457-4, Human PMS1 Recombinant Matched Antibody Pair - PBS only. 83456-8-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag28132. 83456-7-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL



Jurkat cell lysate was measured. The human PMS1 concentration of detected samples was determined to be 6.03 ng/mL (based on a 1.0 mg/mL extract load).