

For Research Use Only

LPAR4 Recombinant antibody

Catalog Number: 83457-1-RR



Basic Information

Catalog Number:	83457-1-RR	GenBank Accession Number:	BC069996	Purification Method:	Protein A purification
Size:	1000 µg/ml	GeneID (NCBI):	2846	CloneNo.:	240079C3
Source:	Rabbit	UNIPROT ID:	Q99677		
Isotype:	IgG	Full Name:	lysophosphatidic acid receptor 4		
Immunogen Catalog Number:	AG17817	Calculated MW:	370 aa, 42 kDa		

Applications

Tested Applications:
ELISA
Species Specificity:
Human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 83457-1-RR against Human LPAR4 were performed. The affinity constant is below 1 pM.