

For Research Use Only

GRIN2A Recombinant antibody, PBS Only (Detector)

Catalog Number: 83465-1-PBS



Basic Information

Catalog Number:

83465-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG29101

GenBank Accession Number:

NM_000833

GeneID (NCBI):

2903

UNIPROT ID:

Q12879

Full Name:

glutamate receptor, ionotropic, N-methyl D-aspartate 2A

Calculated MW:

165 kDa

Purification Method:

Protein A purification

CloneNo.:

240445B5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00461-3, GRIN2A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83465-4-PBS. Detection antibody: 83465-1-PBS. Standard: Ag29101. Range: 0.313-40 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83465-1-PBS against Human GRIN2A was performed. The affinity constant is 21.9 nM.