

For Research Use Only

# HSPA2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83466-2-PBS

## Basic Information

Catalog Number: 83466-2-PBS	GenBank Accession Number: BC036107	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3306	CloneNo.: 240421B4
Source: Rabbit	UNIPROT ID: P54652	
Isotype: IgG	Full Name: heat shock 70kDa protein 2	
Immunogen Catalog Number: AG3539	Calculated MW: 639 aa, 70 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00465-2, HSPA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83466-1-PBS. Detection antibody: 83466-2-PBS. Standard: Ag3539. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00465-3, HSPA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83466-3-PBS. Detection antibody: 83466-2-PBS. Standard: Ag3539. Range: 1.25-80 ng/mL