

For Research Use Only

# GSTT2B Recombinant antibody, PBS Only (Detector)



Catalog Number: 83469-1-PBS

## Basic Information

Catalog Number: 83469-1-PBS	GenBank Accession Number: BC071700	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 653689	CloneNo.: 240493D2
Source: Rabbit	UNIPROT ID: POCG30	
Isotype: IgG	Full Name: glutathione S-transferase theta 2B (gene/pseudogene)	
Immunogen Catalog Number: AG11828	Calculated MW: 244 aa, 28 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00473-1, GSTT2B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83469-2-PBS. Detection antibody: 83469-1-PBS. Standard: Ag11828. Range: 0.78-100 ng/mL