

For Research Use Only

# GSTT2B Recombinant antibody, PBS Only



Catalog Number: 83469-4-PBS

## Basic Information

Catalog Number: 83469-4-PBS	GenBank Accession Number: BC071700	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 653689	CloneNo.: 240493G4
Source: Rabbit	UNIPROT ID: P0CG30	
Isotype: IgG	Full Name: glutathione S-transferase theta 2B (gene/pseudogene)	
Immunogen Catalog Number: AG11828	Calculated MW: 244 aa, 28 kDa	

## Applications

Tested Applications:  
WB, ELISA

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83469-4-RR (GSTT2B antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83469-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83469-4-RR against Human GSTT2B were performed. The affinity constant is below 1 pM.