

For Research Use Only

# IKZF4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83471-1-PBS

## Basic Information

Catalog Number:

83471-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34767

GenBank Accession Number:

BC160169

GeneID (NCBI):

64375

UNIPROT ID:

Q9H2S9

Full Name:

IKAROS family zinc finger 4 (Eos)

Purification Method:

Protein A purification

CloneNo.:

240357C8

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00475-2, IKZF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83471-4-PBS. Detection antibody: 83471-1-PBS. Standard: Ag34767. Range: 0.625-40 ng/mL



Cytometric bead array standard curve of MP00475-3, IKZF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83471-2-PBS. Detection antibody: 83471-1-PBS. Standard: Ag34767. Range: 0.625-40 ng/mL