

For Research Use Only

IKZF4 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83471-2-PBS

Basic Information

Catalog Number: 83471-2-PBS	GenBank Accession Number: BC160169	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 64375	CloneNo.: 240357D3
Source: Rabbit	UNIPROT ID: Q9H2S9	
Isotype: IgG	Full Name: IKAROS family zinc finger 4 (Eos)	
Immunogen Catalog Number: AG34767		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00475-1, IKZF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83471-2-PBS. Detection antibody: 83471-3-PBS. Standard: Ag34767. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00475-3, IKZF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83471-2-PBS. Detection antibody: 83471-1-PBS. Standard: Ag34767. Range: 0.625-40 ng/mL