

For Research Use Only

MIG6; ERFFI1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83475-3-PBS



Basic Information

Catalog Number:

83475-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33341

GenBank Accession Number:

NM_018948

GeneID (NCBI):

54206

UNIPROT ID:

Q9UJM3

Full Name:

ERBB receptor feedback inhibitor 1

Calculated MW:

51kd

Purification Method:

Protein A purification

CloneNo.:

240447A4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00464-2, MIG6; ERRFI 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83475-1-PBS. Detection antibody: 83475-3-PBS. Standard: Ag33341. Range: 0.625-80 ng/mL.