

For Research Use Only

# FOXC2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83476-3-PBS



## Basic Information

Catalog Number:

83476-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33121

GenBank Accession Number:

NM\_005251

GeneID (NCBI):

2303

UNIPROT ID:

Q99958

Full Name:

forkhead box C2 (MFH-1,  
mesenchyme forkhead 1)

Calculated MW:

54kd

Purification Method:

Protein A purification

CloneNo.:

240403C6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00470-1, FOXC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83476-3-PBS. Detection antibody: 83476-1-PBS. Standard: Ag33121. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00470-2, FOXC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83476-3-PBS. Detection antibody: 83476-4-PBS. Standard: Ag33121. Range: 0.313-40 ng/mL