

For Research Use Only

# Mouse MMR/CD206 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83485-4-PBS



## Basic Information

Catalog Number:

83485-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_008625.2

GeneID (NCBI):

17533

UNIPROT ID:

Q2HZ94

Full Name:

mannose receptor, C type 1

Calculated MW:

165 kDa

Purification Method:

Protein A purification

CloneNo.:

240642C8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00623-1, MOUSE CD206 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83485-4-PBS. Detection antibody: 83485-3-PBS. Standard: Eg0790. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP00623-2, MOUSE CD206 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83485-4-PBS. Detection antibody: 83485-2-PBS. Standard: Eg0790. Range: 0.78-100 ng/mL.