

For Research Use Only

FAM38B Recombinant antibody, PBS Only (Detector)

Catalog Number: 83488-2-PBS



Basic Information

Catalog Number:

83488-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24528

GenBank Accession Number:

AB527139

GeneID (NCBI):

63895

UNIPROT ID:

Q9H5I5

Full Name:

family with sequence similarity 38, member B

Calculated MW:

318 kDa

Purification Method:

Protein A purification

CloneNo.:

240406D8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00488-1, FAM38B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83488-3-PBS. Detection antibody: 83488-2-PBS. Standard: Ag24528. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00488-3, FAM38B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83488-4-PBS. Detection antibody: 83488-2-PBS. Standard: Ag24528. Range: 0.313-40 ng/mL