

For Research Use Only

HTRA2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83492-2-PBS

Basic Information

Catalog Number:	83492-2-PBS	GenBank Accession Number:	BC000096	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	27429	CloneNo.:	240496B3
Source:	Rabbit	UNIPROT ID:	O43464		
Isotype:	IgG	Full Name:	HtrA serine peptidase 2		
Immunogen Catalog Number:	AG8455	Calculated MW:	49 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00491-3, HTRA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83492-2-PBS. Detection antibody: 83492-4-PBS. Standard: Ag8455. Range: 0.781-100 ng/mL.