

For Research Use Only

LAMR1/RPSA Recombinant antibody, PBS Only (Capture)

Catalog Number: 83495-3-PBS



Basic Information

Catalog Number:

83495-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6033

GenBank Accession Number:

BC050688

GeneID (NCBI):

3921

UNIPROT ID:

P08865

Full Name:

ribosomal protein SA

Calculated MW:

33 kDa

Purification Method:

Protein A purification

CloneNo.:

240461D2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00494-2, LAMR1/RPSA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83495-3-PBS. Detection antibody: 83495-2-PBS. Standard: Ag6033. Range: 0.313-40 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83495-3-PBS against Human LAMR1/RPSA was performed. The affinity constant is 8.64 nM.