

For Research Use Only

SLC25A48 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83498-2-PBS

Basic Information

Catalog Number: 83498-2-PBS	GenBank Accession Number: XM_017009087	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 153328	CloneNo.: 240155G10
Source: Rabbit	UNIPROT ID: Q6ZT89	
Isotype: IgG	Full Name: Solute carrier family 25 member 48	
Immunogen Catalog Number: AG34514	Calculated MW: 311 aa, 33 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00485-1, SLC25A48 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83498-2-PBS. Detection antibody: 83498-1-PBS. Standard: Ag34514. Range: 0.156-20 ng/mL.