

For Research Use Only

UPK3A Recombinant antibody, PBS Only (Detector)



Catalog Number: 83499-1-PBS

Basic Information

Catalog Number:

83499-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_006953

GeneID (NCBI):

7380

UNIPROT ID:

O75631

Full Name:

uroplakin 3A

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

240488H1

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00484-1, UPK3A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83499-2-PBS. Detection antibody: 83499-1-PBS. Standard: SY01289. Range: 0.195-100 ng/mL.