

For Research Use Only

# Mouse Flt3 ligand Recombinant antibody, PBS Only (Capture)

Catalog Number: 83502-2-PBS



## Basic Information

Catalog Number:

83502-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_013520.3

GeneID (NCBI):

14256

UNIPROT ID:

P49772-1

Full Name:

FMS-like tyrosine kinase 3 ligand

Calculated MW:

26KD

Purification Method:

Protein A purification

CloneNo.:

240414D2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00481-1, MOUSE Flt3 ligand Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83502-2-PBS. Detection antibody: 83502-1-PBS. Standard: Eg0677. Range: 0.156-20 ng/mL