

For Research Use Only

FBXO38 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83509-4-PBS

Basic Information

Catalog Number:	83509-4-PBS	GenBank Accession Number:	XM_005268513.1	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	81545	CloneNo.:	240513G5
Source:	Rabbit	UNIPROT ID:	Q6PIJ6		
Isotype:	IgG	Full Name:	F-box protein 38		
Immunogen Catalog Number:	AG35444	Calculated MW:	134kDa, 1188aa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00502-2, FBXO38 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83509-4-PBS. Detection antibody: 83509-2-PBS. Standard: Ag35444. Range: 0.781-100 ng/mL.