

For Research Use Only

MED21 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83510-2-PBS

Basic Information

Catalog Number: 83510-2-PBS	GenBank Accession Number: BC008380	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9412	CloneNo.: 240419E12
Source: Rabbit	UNIPROT ID: Q13503	
Isotype: IgG	Full Name: mediator complex subunit 21	
Immunogen Catalog Number: AG8568	Calculated MW: 144 aa, 16 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00506-1, MED21 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83510-1-PBS. Detection antibody: 83510-2-PBS. Standard: Ag8568. Range: 1.25-40 ng/mL.